# GenoType MTBDR*plus*

VER 2.0

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**IFU-304A-06**

**CE**

**IVD**

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### Diagnostic performance

#### Pulmonary clinical specimens

Diagnostic performance characteristics of the **GenoType MTBDRplus** VER 2.0 were determined in a study [12] with 338 specimens (including sputum, bronchoalveolar lavages, and pleural aspirates) compared to culture (successful cultivation on Loewenstein-Jensen solid medium or in MGIT (BD Diagnostics, Franklin Lakes, USA)) and subsequent speciation using the **GenoType Mycobacterium CM** VER 1.0) and phenotypic drug susceptibility testing (DST). Additionally, the samples were examined by microscopy. Clinical data of the patients were included in the evaluation.

The study site was located in a high MDR-TB burden country. Microscopy and cultivation methods were conducted on site. Aliquots of the NALC-decontaminated sputum specimens were shipped to a second laboratory to perform DNA extraction and the **GenoType MTBDRplus**. Manual DNA extraction was performed with the **GenoLyse**® kit (162 of the 338 sputum samples), automated DNA extraction was carried out on the **GenoExtract**® instrument using the **GXT DNA/RNA Extraction Kit** (176 of the 338 sputum samples) according to the respective instructions for use.

A congruent **GenoType MTBDRplus** positive, clinical positive result was defined either by positivity of culture and **GenoType MTBDRplus** or when only **GenoType MTBDRplus** was positive and culture negative, but TB was indicated by previous culture-based findings of the respective patient. A discrepant result (**GenoType MTBDRplus** positive and culture negative) does not exclude in all cases a TB infection of the patient as for some patients' histories of a probable TB infection were not available.

**Table 1:** Performance characteristics of the **GenoType MTBDRplus** for detection of MTBC from pulmonary clinical specimens compared to culture/**GenoType Mycobacterium CM** (GT Myco CM) and clinical findings

		Smear-positive		Sens: 100% Spec: /* PPV: 100% NPV: /*	Smear-negative		Sens: 80.3% Spec: 98.4% PPV: 98.0% NPV: 83.3%		
		Culture/GT Myco CM and clinic			Culture/GT Myco CM and clinic				
		Positive	Negative		Positive	Negative			
<b>GenoLyse</b> ®	<b>GenoType MTBDRplus</b>	Positive	39	0	Sens: 97.5% Spec: /* PPV: 97.5% NPV: /*	Positive	49	1	Sens: 78.3% Spec: 96.0% PPV: 94.0% NPV: 84.7%
		Negative	0	1		Negative	12	60	
<b>GXT</b>	<b>GenoType MTBDRplus</b>	Positive	39	1	Sens: 96.3% Spec: 100% PPV: 100% NPV: 92.3%	Positive	47	3	Sens: 86.2% Spec: 100% PPV: 100% NPV: 71.4%
		Negative	1	0		Negative	13	72	

Sens, diagnostic sensitivity; Spec, diagnostic specificity; PPV, positive predictive value; NPV, negative predictive value

\* no value due to low sample number

For evaluation of resistance detection, the 156 samples (78 **GenoLyse**® isolates and 78 **GXT** isolates) were used which were MTBC-positive in both culture and **GenoType MTBDRplus**.

**Table 2:** Performance characteristics of the **GenoType MTBDRplus** for detection of RMP resistance from pulmonary clinical specimens compared to culture/DST

		Smear-positive		Sens: 100% Spec: 92.3% PPV: 96.2% NPV: 100%	Smear-negative		Sens: 96.0% Spec: 93.3% PPV: 96.0% NPV: 93.3%		
		Culture/DST			Culture/DST				
		RMP-R	RMP-S		RMP-R	RMP-S			
<b>GenoLyse</b> ®	<b>GenoType MTBDRplus</b>	RMP-R	25	1	Sens: 96.3% Spec: 100% PPV: 100% NPV: 92.3%	RMP-R	24	1	Sens: 86.2% Spec: 100% PPV: 100% NPV: 71.4%
		RMP-S	0	12		RMP-S	1	14	
<b>GXT</b>	<b>GenoType MTBDRplus</b>	RMP-R	26	0	Sens: 96.7% Spec: 87.5% PPV: 96.7% NPV: 87.5%	RMP-R	25	0	Sens: 90.6% Spec: 90.0% PPV: 96.6% NPV: 90.0%
		RMP-S	1	12		RMP-S	4	10	

Sens, diagnostic sensitivity; Spec, diagnostic specificity; PPV, positive predictive value; NPV, negative predictive value;

RMP-R, resistant to rifampicin; RMP-S, sensitive to rifampicin

**Table 3:** Performance characteristics of the **GenoType MTBDRplus** for detection of INH resistance from pulmonary clinical specimens compared to culture/DST

		Smear-positive		Sens: 96.7% Spec: 87.5% PPV: 96.7% NPV: 87.5%	Smear-negative		Sens: 96.7% Spec: 90.0% PPV: 96.6% NPV: 90.0%		
		Culture/DST			Culture/DST				
		INH-R	INH-S		INH-R	INH-S			
<b>GenoLyse</b> ®	<b>GenoType MTBDRplus</b>	INH-R	29	1	Sens: 100% Spec: 100% PPV: 100% NPV: 100%	INH-R	29	1	Sens: 90.6% Spec: 90.0% PPV: 93.5% NPV: 62.5%
		INH-S	1	7		INH-S	1	9	
<b>GXT</b>	<b>GenoType MTBDRplus</b>	INH-R	28	0	Sens: 96.7% Spec: 87.5% PPV: 96.7% NPV: 87.5%	INH-R	29	2	Sens: 90.6% Spec: 90.0% PPV: 93.5% NPV: 62.5%
		INH-S	0	11		INH-S	3	5	

Sens, diagnostic sensitivity; Spec, diagnostic specificity; PPV, positive predictive value; NPV, negative predictive value;

INH-R, resistant to isoniazid; INH-S, sensitive to isoniazid

#### Cultured material

The diagnostic performance characteristics of the **GenoType MTBDRplus** were determined in a study with 74 cultured samples compared to **GenoType Mycobacterium CM** VER 1.0 and phenotypic drug susceptibility testing (DST). The study site was located in a low MDR-TB burden country. Manual DNA extraction was performed using the **GenoLyse**<sup>®</sup> kit according to the instructions for use. From 74 cultures, 49 were positive for *M. tuberculosis* complex (MTBC) and 25 cultures showed growth of nontuberculous mycobacteria. Hence, for resistance detection in cultured material, 49 isolates were available.

**Table 4:** Performance characteristics of the **GenoType MTBDRplus** for detection of MTBC from cultured material compared to culture/**GenoType Mycobacterium CM** (GT Myco CM)

	Culture/GT Myco CM		Sens: 100% Spec: 100% PPV: 100% NPV: 100%
	Positive	Negative	
<b>GenoType MTBDRplus</b>	Positive	49	
	Negative	0	

Sens, diagnostic sensitivity; Spec, diagnostic specificity; PPV, positive predictive value; NPV, negative predictive value

**Table 5:** Performance characteristics of the **GenoType MTBDRplus** for detection of RMP resistance from cultured material compared to culture/DST

	Culture/DST		Sens: /* Spec: 100% PPV: /* NPV: 100%
	RMP-R	RMP-S	
<b>GenoType MTBDRplus</b>	RMP-R	0	
	RMP-S	0	

Sens, diagnostic sensitivity; Spec, diagnostic specificity; PPV, positive predictive value; NPV, negative predictive value; RMP-R, resistant to rifampicin; RMP-S, sensitive to rifampicin

\* no value due to low sample number

**Table 6:** Performance characteristics of the **GenoType MTBDRplus** for detection of INH resistance from cultured material compared to culture/DST

	Culture/DST		Sens: /* Spec: 100% PPV: /* NPV: 100%
	INH-R	INH-S	
<b>GenoType MTBDRplus</b>	INH-R	3	
	INH-S	0	

Sens, diagnostic sensitivity; Spec, diagnostic specificity; PPV, positive predictive value; NPV, negative predictive value; INH-R, resistant to isoniazid; INH-S, sensitive to isoniazid

\* no value due to low sample number

#### Analytical performance

##### Analytical specificity

The specificity of the **GenoType MTBDRplus** VER 2.0 is ensured by the accurate design of specific primers and probes which considers, among others, homology comparisons of the sequences published in gene databases, and by stringent reaction conditions.

The analytical specificity was determined with 61 DNA isolates including the following MTBC strains: *M. tuberculosis*, *M. africanum*, *M. bovis*, *M. canettii*, *M. microti*, and *M. pinnipedii* (all RMP- and INH-sensitive). The following strains not detectable with the test system were analyzed: *Actinomyces naeslundii*, *Aggregatibacter actinomycetemcomitans*, *Bacillus cereus*, *Corynebacterium ammoniagenes*, *C. bovis*, *C. durum*, *Escherichia coli*, *Gordona rubropertinctus*, *Klebsiella oxytoca*, *Mycobacterium abscessus*, *M. alvei*, *M. asiaticum*, *M. avium*, *M. celatum*, *M. chelonae*, *M. chimaera*, *M. fortuitum* (2 sequevars), *M. frederiksbergense*, *M. gastri*, *M. genavense*, *M. goodii*, *M. gordonae*, *M. heckeshornense*, *M. immunogenum*, *M. interjectum*, *M. intermedium*, *M. intracellulare*, *M. lentiflavum*, *M. marinum*, *M. mucogenicum*, *M. palustre*, *M. peregrinum*, *M. scrofulaceum*, *M. shimoidei*, *M. simiae*, *M. smegmatis*, *M. szulgai*, *M. triplex*, *M. ulcerans*, *M. xenopi*, MRSA, *Nocardia abscessus*, *N. africana*, *N. amarae*, *N. asteroides*, *N. farcinica*, *Porphyromonas gingivalis*, *Prevotella intermedia*, *Rhodococcus erythropolis*, *Saccharomonospora glauca*, *Tannerella forsythia*, *Treponema denticola*, *Tsukamurella inchonensis*, *T. pulmonis*.

The six MTBC isolates were correctly identified as RMP- and INH-sensitive MTBC strains. All other 55 isolates displayed no TUB band and no evaluable band pattern for RMP and INH resistances. Hence, an analytical specificity of 100% was achieved.

##### Analytical sensitivity

For determination of analytical sensitivity of the **GenoType MTBDRplus** for clinical samples, ten parallel BCG cultures were diluted and spiked into MTBC-negative sputum samples (final concentrations in sputum samples: 1000, 500, 160, and 100 bacteria/ml). Including a negative control, DNA was extracted once using the **GenoLyse**<sup>®</sup> kit and once using the **GXT DNA/RNA Extraction Kit**, and analyzed with the **GenoType MTBDRplus** applying the "MDR DIR" PCR protocol. A limit of detection of 160 bacteria/ml was determined with both extraction methods.

For determination of analytical sensitivity of the **GenoType MTBDRplus** for culture samples, four BCG cultures (RMP- and INH-sensitive,  $1.6 \times 10^4$ ,  $1.6 \times 10^3$ ,  $1.6 \times 10^2$ , and 100 bacteria/ml) were set up in triplicate. Including a negative control, DNA was extracted using the **GenoLyse**<sup>®</sup> kit and analyzed with the **GenoType MTBDRplus** applying the "MDR CUL" PCR protocol. A limit of detection of  $1.6 \times 10^4$  bacteria/ml was determined.

#### Reproducibility

##### Intra-assay precision

In order to determine the intra-assay precision of the **GenoType MTBDRplus**, two BCG cultures (RMP- and INH-sensitive, 1,500 and 150 bacteria/ml) and an *M. avium* culture (10,000 bacteria/ml) were set up in triplicate and spiked into negative sputum specimens. These samples and a negative control were tested with the **GenoType MTBDRplus** under identical conditions, applying the "MDR DIR" PCR protocol. DNA extraction was performed once using the **GenoLyse**<sup>®</sup> DNA extraction kit, and once using the **GXT DNA/RNA Extraction Kit**. All parallels showed identical and correct banding patterns and comparable signal strengths. Additionally, signal strengths between the two DNA extraction methods and between different dilutions of the same samples were comparable. Hence, an intra-assay precision of 100% was achieved.

#### Inter-assay precision

In order to determine the inter-assay precision of the **GenoType MTBDRplus**, two BCG cultures (RMP- and INH-sensitive, 1,500 and 150 bacteria/ml) and an *M. avium* culture (10,000 bacteria/ml) were set up in triplicate and spiked into negative sputum specimens. These samples and a negative control were tested in nine runs: on three different days, using three different sets of instruments, and conducted by three different operators. DNA extraction was performed once using the **GenoLyse**® DNA extraction kit and once using the **GXT DNA/RNA Extraction Kit**. The "MDR DIR" PCR protocol was applied for PCR. Apart from the varied parameter, all other testing conditions were identical. No deviations were detected between parallel samples, that is between runs banding patterns were identical and correct and signal strengths were comparable. Moreover, signal strengths were comparable between different DNA extraction methods and different bacterial concentrations. Hence, the inter-assay precision was 100%.

#### Interfering substances

There are substances that may inhibit PCR reactions. Such inhibitors may, for example, originate from the culture medium. In order to assess if the medium influences the **GenoType MTBDRplus**, 6 different *M. tuberculosis* complex samples (4 RMP- and INH-resistant, 2 RMP- and INH-sensitive) were cultured in 4 different media (solid media: Loewenstein-Jensen, Stonebrink, and Middlebrook-7H10, liquid medium: MGIT (BD Diagnostics, Franklin Lakes, USA)). DNA was extracted from the culture samples using the **GenoLyse**® DNA extraction kit and then tested with the **GenoType MTBDRplus**.

All *M. tuberculosis* complex samples showed the same correct results. Hence, it can be excluded that the tested media import inhibitors into the **GenoType MTBDRplus** test.

Interfering substances may also be carried over from the sample material. Hence, the substances indicated in table 7 were tested in order to assess a potential interference of the **GenoType MTBDRplus**. Defined BCG culture dilutions above and at the detection limit were spiked with various amounts of the potential inhibitors. From all samples, DNA extraction was performed using the **GenoLyse**® DNA extraction kit. Then the culture dilutions were tested with the **GenoType MTBDRplus**.

**Table 7:** Tested potential interferents of the **GenoType MTBDRplus**

Substance/class	Description/active ingredient	Test concentration(s)
Allergy relief medicines	Tea tree oil	0.008 % v/v to 0.5 % v/v
Anesthetics (endotracheal intubation)	Lidocaine HCl 4%	20% v/v; 30% v/v
Anesthetics (oral)	Benzocaine 20%	5% w/v
Antibiotics (nasal ointment)	Mupirocin	1.2 mg/ml; 2.4 mg/ml
Antibiotics (systemic)	Amoxicillin	2.2 µg/ml
Anti-tuberculosis drugs	Isoniazid 1 mg/ml	50 µg/ml
Anti-tuberculosis drugs	Rifampicin 1 mg/ml	25 µg/ml
Anti-tuberculosis drugs	Pyrazinamide 10 mg/ml	100 µg/ml
Anti-tuberculosis drugs	Ethambutol 1 mg/ml	5 µg/ml; 50 µg/ml
Anti-tuberculosis drugs	Streptomycin 1 mg/ml	25 µg/ml
Anti-viral drugs	Zanamivir	800 µg/ml
Blood	Whole blood	0.2% v/v to 5% v/v
Blood	Hemoglobin	0.05% v/v to 0.3% v/v
Bronchodilators	Theophylline	222 pmol/ml
DNA (human)		10 µg/ml
Expectorants (oral)	Guaifenesin 400 mg/pill	2.5 mg/ml; 5 mg/ml
Gastric acid	0.5% HCl, 0.1 M KCl, 0.1 M NaCl, pH 1-2	5% v/v
Influenza vaccine as nasal spray (FluMist®)	Live attenuated influenza vaccine	5% v/v
Inhaled bronchodilators	Salbutamol sulfate 2.5 mg/3 ml	50 µg/ml; 100 µg/ml
Mouthwash/gargle solutions	Listerine (eucalyptol 0.029%, menthol 0.042%, methyl-salicylate 0.06%, thymol 0.064%, denatured alcohol 20%)	20% v/v
Mucin: bovine submaxillary gland, type I-S	Purified mucin protein 5% w/v	1.5% w/v; 5% w/v
Nasal corticosteroids	Dexamethasone	1.52 pmol/ml
Nasal gel (homeopathic)	Sulfur	5% w/v
Nasal sprays or drops	Phenylephrine 0.5%	25% v/v; 100% v/v
Nebulizing solutions (hypertonic saline)	NaCl	3% w/v; 5% w/v
Physiologic saline	NaCl (0.9%)	0.9% w/v
<i>Pneumocystis jiroveci</i> medications	Pentamidine	300 ng/ml
Pus		0.2% v/v to 5% v/v
Specimen processing reagents	NALC-NaOH (N-acetyl-L-cysteine-sodium hydroxide)	0.5% v/v; 1% v/v
Tobacco	Nicogel (40% tobacco extract)	0.5% w/v

Inhibition of the **GenoType MTBDRplus** VER 2.0 (invalid test result) was observed in the presence of the following substances in the concentrations as indicated: whole blood at 0.6%, hemoglobin at 0.1%, pus at 2%, lidocaine at 30%, mupirocin at 2.4 mg/ml, tea tree oil at 0.5%, and guaifenesin at 5 mg/ml.

#### Stability

Shelf life of the **GenoType MTBDRplus** test kit when stored as recommended: see box label.

Stability is determined according to DIN EN ISO 23640.

## References

1. World Health Organization. Global tuberculosis report 2014. WHO/HTM/TB/2014.08. World Health Organization, Geneva, Switzerland 2014.
2. Zhang Y, Yew WW. Mechanisms of drug resistance in *Mycobacterium tuberculosis*. *Int J Tuberc Lung Dis* 2009; 13: 1320–1330.
3. Biosafety in microbiological and biomedical laboratories, 5th edition. U.S. Department of Health and Human Services, Centers for Disease Control and Prevention, Atlanta, USA 2009.
4. Protection of laboratory workers from occupationally acquired infections. Approved guideline. Clinical and Laboratory Standards Institute (formerly National Committee for Clinical Laboratory Standards), USA, Document M29 (please refer to the latest version).
5. Kent PT, Kubica GP. Public health mycobacteriology: a guide for the level III laboratory. U.S. Department of Health and Human Services, Centers for Disease Control and Prevention, Atlanta, USA 1985.
6. Isenberg HD. Clinical microbiology procedures handbook. American Society for Microbiology, Washington, D.C., USA 1992.
7. Richter E, Beer J, Diel R, Hillemann D, Hoffmann H, Klotz M, Mauch H, Rüsche-Gerdes S. MiQ 5, Tuberkulose, Mykobakteriose. In: Podbielski A, Herrmann M, Kniehl E, Mauch H, Rüssmann H (eds): Mikrobiologisch-infektiologische Qualitätsstandards. Elsevier, Munich, Germany 2010.
8. DIN, Deutsches Institut für Normung e.V. (ed). DIN 58943-4:2009-02: Medical microbiology - Diagnosis of tuberculosis - Part 4: Primary samples for the diagnosis of tuberculosis and mycobacteria – Qualitative and quantitative requirements, extraction, transport and storage. Beuth, Berlin, Germany 2009.
9. Telenti A, Imboden P, Marchesi F, Lowrie D, Cole S, Colston MJ, Matter L, Schopfer K, Bodmer T. Detection of rifampicin-resistance mutations in *Mycobacterium tuberculosis*. *Lancet* 1993, 341: 647-650.
10. Alonso M, Palacios JJ, Herranz M, Penedo A, Menéndez A, Bouza E, García de Viedma D. Isolation of *Mycobacterium tuberculosis* strains with a silent mutation in *rpoB* leading to potential misassignment of resistance category. *J Clin Microbiol* 2011; 49: 2688-2690.
11. Musser JM. Antimicrobial agent resistance in mycobacteria: molecular genetic insights. *Clin Microbiol Rev* 1995; 8: 496-514.
12. Crudu V, Stratan E, Romancenco E, Allerheiligen V, Hillemann A, Moraru N. First evaluation of an improved assay for molecular genetic detection of tuberculosis as well as RMP and INH resistances. *J Clin Microbiol* 2012; Epub ahead of print, doi:10.1128/JCM.05903-11.

## Important Changes in IFU-304A-06

Chapter	Change
Precautions for Handling Kit Constituents	Classification of the kit reagents according to GHS criteria [Regulation (EC) No. 1272/2008]

**Hain Lifescience GmbH**

Hardwiesenstraße 1, 72147 Nehren, Germany  
[www.hain-lifescience.de](http://www.hain-lifescience.de), +49 (0) 74 73- 94 51- 0